

QUICK PROTOCOL BIOSCRIPT REVERSE TRANSCRIPTASE



N - HOT

QUICK PROTOCOL

BIOSCRIPT REVERSE TRANSCRIPTASE

BioLinker's Bioscrip Reverse Transcriptase™ is an engineered version of MMLV RTx with more thermal stability and increased specificity by reduced RNase H activity mutated domain [1]. The enzyme can be used to synthesize first-strand cDNA at temperatures up to 55°C, providing higher yields of cDNA, and more fulllength product than other reverse transcriptases. It can generate cDNA from 100 bp to >12 kb.

BIOSCRIPT REVERSE TRANSCRIPTASE COMPONENTS

- BioScript RTx™: 10,000 U/mL (200 U/μl)
- Buffer RT 5X* 1000 μl

*250 mM Tris-HCl (pH 8.3 at room temperature), 375 mM KCl, 15 mM MgCl₂

**RNA inhibitor not included

UNIT DEFINITION

One unit incorporates 1 nmol of dTTP into acid-precipitable material in 10 min at 37°C using poly(A)•oligo(dT)₂₅ as template-primer (3).

STORAGE BUFFER

Storage Buffer: 40 mM Tris-HCl (pH 7.5), 200 mM NaCl, 100 mM DTT, 0.01% (v/v) tween 20, 50% (v/v) glycerol

STORAGE CONDITIONS

Store all components at -20°C (non-frost-free). Thaw Buffer RT 5x at room temperature just prior to use and refreeze immediately.

References 1. Kotewicz, M.L., D'Alessio, J.M., Driftmier, K.M., Blodgett, K.P., and Gerard, G.F. (1985) Gene 35, 249.

QUICK PROTOCOL BIOSCRIPT REVERSE TRANSCRIPTASE

STEP 01: Compound and Primers Mix

Add the following components to a nuclease-free microcentrifuge tube:

- 1 μ l of oligo(dT)20 (50 μ M); or 200–500 ng of oligo(dT)12-18; or 50–250 ng of random primers; or 2 pmol of gene-specific primer
- 10 pg–5 μ g total RNA or 10 pg–500 ng mRNA
- 1 μ l 10 mM dNTP Mix (10 mM each dATP, dGTP, dCTP and dTTP at neutral pH)
- Sterile, distilled water to 13 μ l

STEP 02: Mixture an Incubation

Heat mixture to 65°C for 5 minutes
Incubate on ice for at least 1 minute

STEP 03: Collect Content and Additions

Collect the contents of the tube by brief centrifugation and add:

4 μ l Buffer 5X RT 1 μ l Recombinant RNase Inhibitor.

Note: When using less than 50 ng of starting RNA, the addition of RNase inhibitor is essential.

1 μ l of BIOScript™ RTx (200 units/ μ l)*

*If generating cDNA longer than 5 kb at temperatures above 50°C using a gene-specific primer or oligo(dT)20, the amount of BioScript RT may be raised to to 400 U (2 μ l) to increase yield.

QUICK PROTOCOL BIOSCRIPT REVERSE TRANSCRIPTASE

STEP 04: Mix and spin

Mix by pipetting gently up and down and spin the tube.
If using random primers, incubate tube at 25°C for 5 minutes.

STEP 05: Incubation

Incubate **at 50°C for 30–60 minutes.**

Increase the reaction temperature **to 55°C for gene-specific primer.**

Reaction temperature **may also be increased to 55°C for difficult templates** or templates with high secondary structure.

STEP 06: Inactivation

Inactivate the reaction by heating at 70°C for 15

The cDNA can now be used as a template for amplification in PCR. However, amplification of some PCR targets (those >1 kb) may require the removal of RNA complementary to the cDNA.

To remove RNA complementary to the cDNA, add 1 µl (2 units) of E. coli RNase H and incubate at 37°C for 20 minutes.

QUICK PROTOCOL

BIOSCRIPT REVERSE TRANSCRIPTASE

REACTION PROTOCOL

Recommended Setup Reaction:

Components	Volume
5X Buffer RT [4 μ L
10 mM dNTP Mix	1 μ L
Sense primer (10 μ M)*	0.5 μ L
Antisense primer (10 μ M)*	0.5 μ L
BIOSCRIPT reverse transcriptase (200U/ μ l)	1 μ L
RNAse Inhibitor	1 μ L
RNA template (20pg-1000 ng)	1-12 μ L
Autoclaved, distilled water ** (complete till 20ul)	- μ L
TOTAL VOLUME	To 20 μL

*Primers concentration should be optimized

** DEPC water and ultrapure autoclaved water are recommended



If you need any help, we are always available for assistance.

Contacts and Support

✉ yourfriends@biolinker.tech

🌐 www.biolinker.tech

📷 [@biolinker_tech](https://www.instagram.com/biolinker_tech)

📘 [@biolinker](https://www.facebook.com/biolinker)

🌐 [BioLinker](https://www.linkedin.com/company/BioLinker)

📞 [+55 \(11\) 3039-8362](tel:+551130398362)