

CoA

Biolinker / Av. Prof. Lineu Prestes, 2242 / São Paulo – SP / 05508-000 / Brazil
info@biolinker.tech / www.biolinker.tech
twitter: @biolinker_tech



SARS-CoV-2 positive control OrfN5

Product Reference: SARSCOV2-CPORFN5

RNA synthetic made by in vitro transcription. (NB1)

Revision date: 05.08.2021
Version: 1.1

Product Data Sheet

Description	SARS-CoV-2 positive control ORFN5
Synonym(s)	covid, covid-19, covid19, coronavirus, Severe acute respiratory Coronavirus positive RNA control, SARS Cov-2 gene, ORFN5, RNA positive control
HS Code:	29349934
Pack Size	10ug
MW	584pb
Form	Liquid
Shipped in	Wet ice
Storage temp.	-20°C.
Expression System	In vitro Transcription
Purity	Not determined
Impurities	rNTPs, enzymes denatured, BSA, salts from buffer-HEPES ph. 7.5
GenBank Accession:	

Lot: 21_03

Certificate of analysis evaluated by integrity and concentration using RNA pico 6000 electropherogram and qubit kits fluorimeter.

Technical responsible